

“ENT MICROSCOPE”



MODEL: ENT-MS/23E



Eyepiece 12.5x/ 10x
0~200° Binocular tube



180-300mm zoom
lens



CCD



CCD adopter



360° Assistant



splitter



Mobile adopter



Camera adopter

Address :S.S.Padana Superiore, 207 - 85624 Vimodrone (MI) ITALY

Tel: +39 02 6532489

“TECHNICAL DATA”

Eyepieces

magnification 12.5x, pupil distance adjustable range 55mm~75mm, diopter adjustment range +6D ~ -6D.

Binocular tube

0~200° binocular tube

Magnification

manual 3-step changer, ratio 0.6, 1.0, 1.6, total magnification 6x, 9x, 14x (F 250mm),

illumination

LED cold light, strong light, longer lifetime, brightness adjustable infinitely, illumination diameter of light spot Φ 70mm, illumination intensity \geq 60000lux, field of view Φ 53mm, Φ 32mm, Φ 20mm build-in small-spot, both blue and yellow light filter.

Objectives

F=250mm (180mm-300mm adjustable objectives for optional)

Rotating head for all direction

operate by two handles, free and easy to operate in 0-360°.

Maximum length of arm 1100mm



Q tech Medical is the first manufacturer who produces operation microscopes in Italy. We have begun to produce ENT operation microscopes since 1970s. And from 1991, we began to manufacture other optical-electronical products, such as spotting scope and photoetching machines. And now you can find optical-electronical products made by **Q tech Medical** all over the world.



Address :S.S.Padana Superiore, 207 - 85624 Vimodrone (MI) ITALY

Tel: +39 02 6532489